



PULSE OXIMETER HANDHELD EST-01 PH



Spesifikasi

- Display Size : 3 Inch (240 x 400 dot)
- Weight : Approx 273 g (Main)
- Dimension : Approx 72 x 154 x 26 mm
- Power Input : 100-240 Vac, 50/60Hz
- Battery : L-ion Rechargeable Battery (2600mA)
- Range of SPO2 Display : Up to 100%
- Resolution : 1%
- Accuracy : Adult : 70-100% \pm 3 Digits
Neonate : 70-100% \pm 3 Digits
- Perfusion Index : 0,05 - 20%
- Measurement Range : Up to 250 bpm
- Resolution : 1 bpm
- Accuracy : \pm 3 Digits
- Memory : Save Continuously For 20 Days
(For 10 Second Saving Period)



PT KARYA PRATAMA

Daan Mogot Arcadia Blok F1 No. 8A-8B Jl. Daan Mogot Km 21,
Kec. Batuceper, Kel. Batuceper, Kota Tangerang, Banten 15122
Telepon : (021) 5595 8085, Fax : (021) 5595 8125